

tweakz® For Aligners

4-in-1 Ortho Care Tool & Travel Case

- 1 Remove retentive aligner trays
- 2 Remove and replace rubber bands
- 3 Remove elastics and debris
- 4 Smooth rough spots that cause ulcerations



INTENDED PURPOSE

Tweakz for Aligners is intended to assist with:

- Aligner Removal
- Difficulty Attaching or Removing Rubber Bands to the Attachments
- Remove Food or Debris Stuck between Teeth or around the attachment
- Smooth Sharp attachment or sharp piece of aligner causing ulceration or irritation around the cheek, lip, or soft tissue.

INDICATION

Users of orthodontic aligners.

INTENDED USERS

Adolescents and adults with orthodontic aligners.

CONTRAINDICATIONS

There are no specific contraindications for these devices.



OrthoNu LLC • 125 Half Mile Road,
Suite 200 • Red Bank, NJ, 07701 USA
whatsnu@orthonu.com • orthonu.com



Self-Care Solutions
You got this!

WARNINGS AND PRECAUTIONS

Please call or go to your orthodontist and team, if you should need any assistance in using the Tweakz appliance. They can explain and show you how each tool can be used so you feel comfortable.

DIRECTIONS

Aligner Removal.

1. Open the tweakz 4-in-1 to rotate open the **1** Aligner Remover attachment and slide it around until it comes to a full stop.
2. Gently place the **Aligner Remover** tool toward the back of the aligner by the first molar on either side.
3. Position at an 90 degree angle between edge of aligner margin and gum line.
4. Slowly push downward towards mouth on upper arch and upwards on lower to remove aligner.

Difficulty Attaching or Removing Rubber Bands to the Attachments.

1. Open the tweakz 4-in-1 to rotate open the **2** Rubber Band Applicator attachment and slide it around until it comes to a full stop.
2. Put the elastic into the deepest part of the **Rubber Band Applicator** where it forms a U and then put onto first attachment.
3. Stretch the rubber-band to second attachment and attach by lifting up and around once rubber band is secure.

Food or Debris Stuck Between Teeth or Around the Attachment.

1. Open the tweakz 4-in-1 to rotate open the **3** Dental Pick attachment and slide it around until it comes to a full stop.
2. Gently place the pick onto the surface of tooth and food/debris interface and gently move outward towards the mouth.

Sharp Attachment or Sharp Piece of Aligner Causing Ulceration or Irritation Around Cheek, Lip or Soft Tissue.

1. Open the tweakz 4-in-1 to rotate open the **4** Diamond Dental File and slide it around until it comes to a full stop.
2. Locate the sharp spot as you move the file in a back and forth direction until smooth and no longer problematic.

Cleaning your Tweakz.

Clean the tweakz 4-in1 tool with soap and water or wipe with alcohol to disinfect.

INSTRUCTIONS FOR REUSE

Do not use the device if it is damaged or shows signs of wear and tear.

STORAGE / HANDLING CONDITIONS

No special storage or handling instructions are necessary.

INSTRUCTIONS FOR DISPOSAL

Dispose of this device in accordance with local regulations and guidelines.

Any serious incident that has occurred in relation to the use of the device should be reported to the Manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Date of issue: 04-DEC-2024

Revision: 1.0



MedEnvoy Global B.V.
Prinses Margrietplantsoen
33 – Suite 123, 2595 AM
The Hague, The Netherlands



MedEnvoy Switzerland
Gotthardstrasse 28,
6302 Zug, Switzerland

Australian Sponsor:
MedEnvoy Global Pty Ltd
60 Martin Place
Sydney, NSW 2000
Australia